



Human Serum Albumin (HSA)

For laboratory procedures only; other uses must be qualified by the end user.

| Product Description | REF Number | Unit Size |
|--------------------------------|------------|-----------|
| Human Serum Albumin, 100 mg/mL | ART-3001 | 12 x 5 mL |
| Human Serum Albumin, 100 mg/mL | ART-3003 | 100 mL |

INTENDED USE

Human Serum Albumin (HSA) has been designed to be used as a tissue culture protein supplement in those applications requiring additional supplementation. Protein in the form of albumin, which makes up the bulk of the protein present in blood serum and reproductive tract fluid, is added to tissue culture media because it is thought to maintain the stability of cell membranes and chelate trace amounts of toxic components present in the culture water, media components, and containers and culture dishes. Protein, in the form of patient's serum or albumin, has been used extensively in media for IVF, GIFT, ICSI, embryo culture, embryo transfer, cryopreservation, and sperm washing for IUI.

DESCRIPTION

Product contains 100 mg/mL total protein (weight/volume) in saline solution. Each lot of HSA is tested for pH (7.4 ± 0.2), osmolality (280 ± 10 Osm/kg water), sterility (no detectable contamination), and biocompatibility (>80% mouse zygote development to blastocysts). All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HB_sAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2°C to 8°C. Warm to ambient or incubator (37°C) temperature prior to use. Do not freeze or expose to temperatures greater than 39°C. The product is stable until the expiration date shown on the label or within 30 days of the Date of First Use provided that proper aseptic procedures have been observed by the user:

- Remove desired volume of product using aseptic procedures.
- Once product has been removed from the original container, reseal the container to ensure a tight seal. Write the date the product was first opened on the product label. Do not use product longer than 30 days after opening the container.
- Once removed, do not return any volume of product to the original container.
- Once the product has been opened, store the sealed container at 2°C to 8°C.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin gel clot tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

SUGGESTED DIRECTIONS FOR USE

For sperm preparation and embryo culture: Use at 5 mg/mL. For 10 mL of medium, add 0.5 mL of HSA solution to 9.5 mL of bicarbonate-buffered medium, e.g., Quinn's Advantage® Fertilization (HTF) Medium (REF #1020/1021) or Quinn's Advantage® Cleavage Medium (REF #1026/1027). *NOTE:* for washed sperm samples for IUI, use QUINN'S® Sperm Washing Medium (REF #1005/1006), which already contains 5 mg/mL HSA.

For embryo transfer: Use at 30 mg/mL. For 10 mL of medium, add 3.0 mL of HSA solution to 7.0 mL of Quinn's Advantage® Medium with HEPES (REF #1023/1024).

For embryo cryopreservation: Use at 12 mg/mL. For 10 mL of medium, add 1.2 mL of HSA solution to 8.8 mL of Quinn's Advantage® Medium with HEPES (REF #1023/1024).

For micromanipulation (ICSI and Assisted Hatching): Use at 5 mg/mL. For 10 mL of medium, add 0.5 mL of HSA solution to 9.5 mL of Quinn's Advantage® Medium with HEPES (REF #1023/1024).

Laboratories may establish through appropriate testing that higher or lower concentrations than those suggested above are optimal for specific applications.

NOT INTENDED FOR INJECTION BY USUAL PARENTERAL ROUTES IN HUMANS OR ANIMALS.

Each laboratory should make its own determination of which medium and protocol to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog (REF #80572).

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter.

To avoid problems with contamination, aseptic technique should be used to remove aliquots of product from the container. Do not mix components of one container with another. Discard minimal amounts of excess product remaining in the bottle.

This product contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for

hepatitis B surface antigen (HB_sAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Caution: Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

RELATED PRODUCTS

SAGE Assisted Reproduction Products™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE Support Line at the number below.

QUINN'S® and Quinn's Advantage® are registered trademarks of CooperSurgical, Inc.

Patrick Quinn, PhD, HCLD, SAGE IVF, Inc.

Call the **SAGE SUPPORT LINE** at:
1-800-243-2974 or 1-203-601-5200



SAGE In-Vitro Fertilization, Inc. A CooperSurgical Company
95 Corporate Drive, Trumbull, CT 06611
www.coopersurgical.com