



SAGE In Vitro Maturation Media

For In Vitro maturation of immature oocytes

Benefits

- Immature oocytes can be cultured to maturity in vitro
- A treatment option for poor responder or over-responders to ovarian stimulation
- Avoids side effects resulting from gonadotropin stimulation including the risk of OHSS
- IVF clinics can offer new treatment options to patients
 - Reduces patient costs
 - Less patient discomfort
 - Fewer drugs



New Treatment Options For Your Patients

Product Overview

In Vitro Maturation is a process by which immature oocytes are harvested from the ovaries and cultured to fully developed embryos suitable for transfer. The SAGE In Vitro Maturation Media, our three-medium product, takes the immature oocytes and can cultivate them to robust embryos eliminating the patient's need for stimulating hormones thus creating a "natural cycle" treatment option.

Product Description

IVM Kit, ART-1600, for the in vitro maturation of human oocytes

- 1 X 50 mL Oocyte Washing Medium
- 1 X 20 mL Oocyte Maturation Medium
- 1 X 5 mL Embryo Maintenance Medium

Clinical Results*

- Cycles: 550
- Oocyte maturation: 75% - 80%
- Clinical pregnancy rate: 35% - 40% per embryo transfer
- Live births: 115



Patient Profile

- Women with PolyCystic Ovarian Syndrome (PCOS)
- Women who want to undergo a "natural cycle" IVF procedure by avoiding stimulating hormones
- Women who want to avoid the potential risk of fertility drug use, such as OHSS patients and patients having undergone repeated stimulation cycles
- Fertility preservation for women undergoing radiation treatment for cancer and who do not have the time for natural stimulation protocol

Chian, R., Lim, J., & Tan, S. (2004, June). State of the Art in In-Vitro Oocyte Maturation. *Current Opinion in Obstetrics & Gynecology*, 16(3), 211-219.

* Data on file at CooperSurgical, Inc.

CooperSurgical

www.coopersurgical.com

To find out more about SAGE In Vitro maturation Media,
or to place an order call 800.243.2974.

The marks bearing the symbol "®" are registered trademarks of CooperSurgical.

SAGE® In Vitro Fertilization
A CooperSurgical Company

95 Corporate Drive, Trumbull, CT 06611
203.601.5200 • 800.243.2974 • fx 800.262.0150
International: 203.601.9818 • fx 203.601.4747

Form Number: 80798 Rev 03/08
Copyright © 2008 CooperSurgical